


# Preliminary Results for Brake Wear Particle Emission PN and PM Concentration and Size Distribution Measurements

M. Moisio, O. Vainio Dekati Ltd., A. Sin, S. Ansari  
ITT, J. Vanhanen, Airmodus Ltd.



50<sup>th</sup> PMP meeting 3-4 April 2019 Brussels

AIRMODUS





- Sales, manufacturing and development of fine particle (aerosols) measurement solutions (automotive/power plant emissions, pharmaceutical, air quality, nanotechnology, defence, science/R&D etc.)
- Core competence: Real-time fine particle size & concentration measurement and sample conditioning technologies
- Privately owned technology spin-off company from TUT Aerosol Physics Lab (1994)
- Located in Kangasala Finland
- Exports ~ 95 % of sales, distributors in ~35 countries worldwide
- ~20 high-know-how employees

# Airmodus Ltd

It's the small things that count

Founded in 2010

Spin-off from the University of Helsinki by a group of aerosol experts specialized in nanoparticle detection

Offers particle detectors capable of detecting nanoparticles starting from 1 nm is diameter

Located in Helsinki, Finland



AIRMODUS

Gas Phase   Nano particles   Particles emitted by Traffic





Pollen

# Goal of the study



- Part 1 of longer measurement campaign, part 2 not yet completed
  - Ford WLTP test cycle
- Determination of brake event repeatability
  - Size distribution, real-time behaviour
- Testing of brake dyno, WLTP cycle and temperature dynamics
- Instrument development
- Determination of system operation and instrument operation
  - ELPI+, CPC a20, eFilter (DC+mass)
  - NOT a full instrument comparison!

# ELPI<sup>®</sup>+: Electrical Low Pressure Impactor

- Number size distribution and concentration
  - Real-time, 10 Hz
- 6 nm - 10  $\mu\text{m}$ 
  - 14 size fractions
  - 100 or 500 with High Resolution ELPI<sup>®</sup>+
- Particles are collected
  - Possibility for chemical analysis on the collected samples
- Wide dynamic range
  - From outdoor air to power plant stack concentrations
- High Temperature (HT-ELPI<sup>®</sup>+)
  - For undiluted sampling a up to 180°C



- Standard gravimetric PM filter holder
- Real-time measurement of PM accumulation on the filter
- Both standard and gravimetric measurements integrated in one instrument
- Fully automated operation



Niemelä et al. 26th CRC Real World Emissions Workshop  
Data courtesy by Ford Motor Company





# Airmodus A20 butanol Condensation Particle Counter

Robust CPC for long term measurements


## SPECS:

One flow design with 50% cut-off at **10 nm** (can be delivered with cut-off from 5-10 nm)

Total aerosol particle number concentration with single particle detection from **1-30 000 #/cc**


Every CPC is calibrated for size cut-off and different particle concentrations up to **100 000 #/cc**. Calibration curve fit provided with the instrument.

Designed to be used as stand-alone, with the Airmodus A10 PSM, or within an integrated measurement system such as DMPS, SMPS or PMP.



# Measurement setup

Filter system: 1° G4, 2° F9, 3° H13




- ITT Dyno operates in light underpressure to avoid PM dust dispersion.
- Two independent ventilation engines
- No temperature and RH control (at the moment)



# Measurement Setup




# Brake event determination




# ELPI+ vs CPC concentration & brake events

## Example with eFilter current




# ELPI+ vs CPC concentration & brake events

## Example with eFilter current and brake temperature






# Repeatability (eFilter)



# ELPI+ vs CPC concentration & brake events




## Size distributions






# Size distributions

- Brake event peaks repeatable and shape rather constant
- Wide size distribution
- Nucleation, accumulation and coarse modes visible
- Mass dominated by  $>1 \mu\text{m}$
- Even number concentration show coarse fraction






# ELPI+ vs. CPC





# CPC vs. eFilter



# Total mass, eFilter collection vs. ELPI<sup>®</sup>+

- eFilter collected 2.2 mg over the whole measurement period (3 WLTP cycles plus some background periods)
  - Flow 10.1 lpm, sampling time ~ 28h 57min(10:18 (8/1/2019) - 15:15 (9/1/2019) )
  - Average PM10 mass, by total sampling time: **0.130 mg/m<sup>3</sup>**
  - Average PM10 mass, divided by actual cycles(~18h): **0.202mg/m<sup>3</sup>**
- ELPI<sup>®</sup>+ average mass concentrations (mg/m<sup>3</sup>):

Cycle no.	PM10	PM2.5	PM1
1	0.230	0.054	0.012
2	0.236	0.055	0.017
3	0.278	0.055	0.015

# Preliminary tests -observations


- Good dyno and temperature operation and results
- Particle instruments working well
- Concentrations in operational range (low for the PM mass)
- ELPI+ vs. CPC PN correlation very good
- ELPI+ vs. eFilter SA/LDSA (current) very good
- ELPI+ distributions very good
- eFilter braking spike resolution good
  
- Dyno system under modifications
- Future work
  - Loss modelling
  - PM10, PM2.5, impactor mass analysis
  - Cycle emission estimations

Thank you!

mikko.moisio@dekati.fi



AIRMODUS



@DekatiLtd



www.dekati.com